


New!
 Ready-to-Use
 HLA Gels



- multichannel compatible
- pre-stained less than 7 min.*
- up to 96 samples/run
- easy to handle reliable


Elchrom Scientific
 Excellence in Electrophoresis

www.elchrom.com

ISO 9001:2008 SWISS **TS**

Product information

The new format of the 96 wells ready-to-use HLA gels consists of 3 gels (1 presentation) mounted on a plastic backing for easy handling. Each gel is provided with 4 x 8 wells plus 4 x 1 well for the Marker.

Using HLA Ready-to-Use Gels is multi-beneficial:

- Easy handling! Gels are 8-channel pipette compatible according to SBS standards.
- Pre-stained! HLA Ready-to-Use Gels are available with GelRed™ or EtBr.
- Save time! With HLA Ready-to-Use Gels up to 96 samples can be run in as little as 7 min.*
- Quality controlled reliability! HLA testing with ready-to-use gels is highly reproducible to avoid time consuming, expensive failures and re-testing.
- Increased flexibility! Cost effective use of gels. Just use what you need.

* using ORIGINS by Elchrom™

Ordering

P/N	Product	Quantity
8110	HLA Ready-to-Use Gel (3x32) EtBr	1 box of 9 presentations
8120	HLA Ready-to-Use Gel (3x32) GelRed™	1 box of 9 presentations
8810	HLA Gel Catamaran (3x32)	1
8820	HLA Gel Catamaran (1x32)	1
2100E	ORIGINS by Elchrom™ (230V, EU-plug)	1
2029E	PSE Electrophoresis Power Supply, 200 V/2000 mA	1

For Research Use Only.
GelRed™ is sold under license from Biotium, Inc. Biotium has the following US patents and foreign equivalents for GelRed™ : US 7, 601, 498, PCT/US2008/013698

Contact

Elchrom Scientific
 Cham, Switzerland
 service@elchrom.com
 Phone +41 41 747 25 50
 Fax +41 41 743 25 36

For detailed information visit
www.elchrom.com